



C-EY-02
Cabinet With Overbed Table

Top table and
front covers
made of HPL

Partition for stuff storage

Metal main frame

Thickness of HPL 6mm

Risers made of aluminum material
on the edges of the dining table

Shock absorber controlled
height adjustable dining table

360° swivel cross lock castors

General & Material Features

Telescopic drawer rails with brake


Aluminium columns
supporting frame connections

Electrostatic powder
painted metal parts

Easy to clean and removable
ABS plastic drawer containers

Monoblock aluminum
foldable dining table

Drawer separator



General & Material Features



Dining Table Tilt Movement



Dining Table Height Movement

Width	Lenght	Height	Bed Table Width	Bed Table Lenght	Bed Table Height Min.	Bed Table Height Max.
59 cm	44,5 cm	86 cm	40 cm	58 cm	84 cm	112 cm



Medikal 2000

Technology for Health

B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr