


2
Linear
actuators



DC-02
Blood Donor Chair

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation



Hand Control Unit

Control Units (Standard)
Battery (Optional)



Foot pressing(Optional)

Control units pictures are representative.




70° Backrest Movement



80° Legrest Movement

Width (Over all)	Lenght (Over all)	Height	Backrest Surface	Sitting Surface	Backrest Angle	Legrest Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
95 cm	185 cm	72 cm	60x80 cm	60x50 cm	70°	80°	230 kg	220 – 240 V AC	2





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr