

LFT-01

PATIENT TRANSPORT LIFT


1
Linear
actuators





Actuators With Min IP
X4 Protection For Quiet
Operation

Ø 120 mm and 65mm
360° Swivel Castors



Hand Control




Control Box
Built-in Battery




Height Adjustable Seating Surface



Pedal to turn feet left and
right



General Features



Height (Over all) A	Width (Over all) B	Height Min. C	Height Max. D	Length Min. E	Length Max F	Safety Working Capacity	Operating Voltage	Number of Actuators
150 cm	130 cm	95 cm	187	71	96,5	150 kg	220–240 V AC	1



B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr