



MA - 02
Drug and Anesthesia Trolley

Colour
Options





Disposal bin with lid
made of ABS material

Seal system for
central locking drawers



Telescopic rails for
easy opening and
closing drawers



Medication boxes
attached on top



Nurse writing table


Optional
Key Locking
Feature



Width	Lenght	Height
82.5 cm	58.5 cm	157.5 cm

5 drawers with following dimensions:

2 pcs 445x545x70mm,
2 pcs 445x545x160mm,
1 pc 445x545x250mm



Side panels and top table
made of polypropylen material

Main frame made of
durable aluminum material

360 ° rotatable
2 castors with brakes
and 2 non-brake wheels

Material Features



B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr