




MN - 03
Mechanical Hospital Bed




Impact Protection Bumpers




Ø125 mm 360 ° Swivel
Cross Lock Castors




Height Adjustable Iv Pole




Easy foldable Side Rails



40-73cm Height Movement




30° Legrest Movement



70° Backrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Number of Crank
109 cm	216 cm	40 cm	73	90x190 cm	90x190x12 cm 28DNS	70°	30°	-	-	230 kg	3





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr