

3
Linear
actuators



PC-01
Podiatry Chair



Actuators With Min
IPX4 Protection
For Quiet Operation

Control Box
Built-in Battery




Foot Control




Legrest Support






70° Backrest Movement




30° Trendelenburg



57-80 cm Height Movement



Calf Movement



90° Top Chassis With
Right And Left Turning

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Backrest Lying Surface	Sitting Surface	Backrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
85 cm	168 cm	57 cm	80 cm	58x85 cm	66x51 cm	70°	30°	-	230 kg	220 – 240 V AC	3



ABS Plastic
 Bottom Chassis
 Cover

Nano ceramic coated
 electrostatic powder
 painted metal parts

Two-pieces
 upholstery with
 polyurethan filling

Ability to fixing
 on the floor

Bottom chassis profiles:
 30x30x2mm
 30x50x2mm

Mechanism profiles:
 30x40x2mm

Upper chassis profiles:
 20x20x2mm

Material Features



B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr