




3
Lifting
columns

+

2
Linear
actuators


PLUS T

ELECTROMECHANIC ICU
HOSPITAL BED WITH THREE COLUMNS



Actuators With Min IPX4
Protection For Quiet
Operation

Accessory Pocket



Hand Control Unit
Nurse Panel
With Function Locking Feature



Control Box
Built-in Battery



Manual CPR Arm

Mechanical Angle Indicator
on Siderails
(backrest-trend angle)



Height Adjustable IV Pole



Ø 150 mm 360 ° Swivel
Central Brake Castors
and ABS Plastic Bottom Chassis

Control units pictures are representative.



File Holder

Oxygen Cylinder Holder



Fifth Caster




Actuators With Min IPX4
Protection For Quiet
Operation



Leg Extension Part

Device Stand





Actuators With Min IPX4
Protection For Quiet
Operation




HPL Lying Surface



Optional Features



Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Lateral Tilt Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	58,5cm	95,5cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	12 °	12 °	21 °	230 kg	220 - 240 VAC	5





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.No:19/A1
Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr