




2
Linear
actuators



PLUS A2

Electromechanic Hospital Bed

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø 125 mm 360 ° swivel
cross lock castors
ABS Plastic Bottom Chassis



Hand Control Unit



Height Adjustable Iv Pole




Manual CPR Arm



Mechanical Angle Indicator
on Siderails (backrest- trend angle)




Control units pictures are representative.




Siderails Control Panel

Sticker options
in different colours
for Head/Footboards
and side rails



X-Ray Cassete Holder


Underbed Light




HPL Lying Surface

Leg Extension Part






70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	42,5	85x190 cm	85x190x12 cm 28DNS	70°	30°	-	-	230 kg	220 – 240 V AC	2





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr