



YM-07
Aluminium Overbed Table



Width	Lenght	Height Min.	Height Max.	Bed Table Width	Bed Table Lenght
44,5 cm	83,5 cm	81,5 cm	116 cm	40 cm	73 cm



Medikal 2000

Technology for Health

B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr